

STATIONARY CONCENTRATOR COMPARISON



MANUFACTURER

React Health

React Health

Caire

Philips

DeVilbiss

Caire

5 LITER
(Side by side)



5L CONCENTRATORS

MODEL		Platinum 5NXG*	Stratus 5	Companion™ *	EverFlo/ EverFlo Q	DeVilbiss 5L 525	AirSep VisionAire™ 5*
SPECIFICATIONS	Continuous Oxygen Flow (max LPM)	5	5	5	5	5	5
	Oxygen Concentration	87-96%	88%-98%	87%-95.5%	90%-96%	87%-96%	87%-95.5%
	Oxygen Outlet Pressure (max psi)	5 PSI	6 PSI - 10 PSI	5 PSI	5.5 PSI	8.5 PSI	9 PSI
	Operating Altitude (feet above sea level)	0 to 13,000	N/A	-1,253 to 9,879	0 to 7,500	N/A	-1250to10,000
SIZE	Size Height x Width x Depth (inches)	24" x 16" x 12"	15.3" x 9" x 23.6"	21.5" x 12.5" x 13.5"	23" x 15" x 9.5"	24.5" x 13.5" x 12"	20.8" x 14.1" x 11.5"
	Weight (pounds)	35	34.1	36	31	36	30
NOISE LEVEL (dBA)		<44	≤45	50	45	48	45
HOMEFILL COMPATIBLE		Yes	No	Unknown	No	Unknown	Yes
WARRANTY		3 years	3 years	3 years	3 years	3 years	3 years

*Companion is a Trademark of Caire, Inc., VisionAire 5 is a Trademark of Caire, Inc.

*The Platinum 5NXG is not currently available. Final specifications to be determined.