

Introducing

PHOENIX™ 5L

The 5 LPM Stationary Oxygen Concentrator



Phoenix 5L

Designed with the patient in mind.

Compact, lightweight, and quiet design for discreet, non-intrusive therapy

Low power consumption helps reduce costs to support budget-conscious care

Integrated internal fire break valve for added safety



PHOENIX™ 5L

Designed with your business in mind.

Reliable performance with patented air separation control system.

Low running temperatures reduce wear on valves and sieve beds to extend life.

Intelligent design with air intake from high on the side and exhaust out the bottom helps prevent recirculation of purged gas.

Easy access to flowmeter and resettable circuit breaker.

SERVICEABILITY

Easy access to internal components by removing just 5 screws— unit remains upright during service.

Operable without the shroud to simplify troubleshooting.

Comprehensive troubleshooting guide designed to minimize repair time and ensure only required parts are replaced.



PHOENIX 5 ALARMS

Power Supply Failure
Start-Up Failure
Low and High Pressure
Low Purity
High Temperature
Low Flow
Power Loss



TYPICAL PARAMETERS

Flow rate (L/min)	0.5 - 5.0
Oxygen Concentration	93% +/- 3%
Rated power	328W
Voltage (V/Hz)	AC 120V / 60
Sound pressure level	<49.5 dBA @ 3LPM
Net weight	35 lbs
Dimensions	20 x 13.7 x 11 inches

MODEL NUMBER

SOC5100 - Phoenix™ 5 Stationary Oxygen Concentrator

REPLACEMENT PART NUMBERS

PRTSOC015 -	Cabinet Filter
PRTSOC016 -	Inlet Filter
PRTSOC007 -	Firebreak Valve
PRTSOC022 -	Filter Grid
PRTSOC014 -	Humidifier Bottle
PRTSOC037 -	Humidifier Bottle Adapter
PRTSOC038 -	Humidifier Bottle Tubing
O2S2007 -	Standard 7' Oxygen Cannula
O2U2007 -	Ultra Soft 7' Oxygen Cannula

ReactHealth.com

© 2025 React Health, all rights reserved.
Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.

MI.MK.0406, Rev A