



Date: April 5, 2023

Travel Compliance Letter for React Health CPAP Devices

To whom it may concern:

React Health’s continuous positive airway pressure (CPAP) machines are medical devices prescribed for use on or by the order of a physician and are cleared for sale by the U.S. Food and Drug Administration, a branch of the Federal Government of the United States of America.

The U.S. Department of Transportation (DOT) Final Rule, “Non-discrimination on the Basis of Disability in Air Travel” (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009, provides important new requirements for the accommodation of passengers with respiratory assistive devices.

Specifically, CPAP devices may be used on board an aircraft, without further testing by the carrier, provided they have been tested for electromagnetic compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M. CPAP devices described in this letter are considered FAA compliant.

React Health Part Number	React Health Device Name
LG2000	Luna II CPAP Device
LG2A00	Luna II APAP Device without standardized Wi-Fi
LG2A01	Luna II APAP Device with standardized Wi-Fi Module
LG2A02	Luna II APAP Device with Cellular Modem
LQX2000	Luna II QX CPAP
LQX2A00	Luna II QX Auto PAP
LG3500	Luna G3 CPAP Device with Cellular Modem – HT
LG3500W	Luna G3 CPAP Device – HT (w/o connectivity)



REACTHEALTH

React Health Part Number	React Health Device Name
LG3600	Luna G3 APAP Device with Cellular Modem – HT
LG36001	Luna G3 APAP Device with standardized Wi-Fi Module- HT
LG36002	Luna G3 APAP Device with Cellular Modem – without HT
LG3600W	Luna G3 APAP Device – HT
LG3700	Luna G3 BPAP 25A Device with Cellular Modem – HT
LG37001	Luna G3 BPAP 25A Device with standardized Wi-Fi Module – HT
LG3700W	Luna G3 BPAP 25A Device – HT
LG3800	Luna G3 BPAP BiLevel 30 VT S/T Device with Cellular Modem – HT
LG38001	Luna G3 BPAP BiLevel 30 VT S/T Device with standardized Wi-Fi Module – HT
LG3800W	Luna G3 Luna G3 BPAP BiLevel 30 VT S/T Device – HT

Please contact React Health Customer Service (863) 226-6285 with any questions.

Best regards,

Chuck Underwood

Chuck Underwood
Director, Quality and Regulatory Affairs
React Health