



React Health's best selling full face mask offers the following features for patient comfort and compliance:

- No forehead brace allows for increased line of sight - even with your glasses on
- Self-adjusting ergonomic cushion design allows for superior seal
- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- 360-degree dual-swivel elbow for active sleepers and allows for sleep in any position
- Air vent diffuser for a quiet night sleep
- Anti-phyxia valve safety feature opens for fresh air in case of machine shut off



REACTHEALTH

Siesta
Full Face Mask



The Siesta Full Face Mask is labeled on the back of the mask.

Allows for ease of identification for troubleshooting and resupply



Mask Part Numbers

SFF1000 - Siesta Full Face Mask - All Sizes Kit

SFF1001 - Siesta Full Face Mask - Small

SFF1002 - Siesta Full Face Mask - Medium

SFF1003 - Siesta Full Face Mask - Large

Replacement Supplies Part Numbers

SFF3001 - Replacement Cushion - Small

SFF3002 - Replacement Cushion - Medium

SFF3003 - Replacement Cushion - Large

SFF2001 - Replacement Headgear



REACTHEALTH

www.reacthealth.com

info@reacthealth.com

Phone: 863-226-6285

Fax: 863-226-6284