



## INTRODUCING Stratus 5

### The 5 LPM Stationary Oxygen Concentrator from React Health

- Slim line design for a smaller footprint in the home
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet operation - compressor is specially placed, and padded to reduce noise
- Simple, intuitive operation



REACTHEALTH



**Stratus 5L**  
Designed with the Patient in Mind

# STRATUS 5L

Designed with your business in mind



## Simple to service

Easy access to internal parts

## 3-year standard warranty

(2 year extended warranty available)

## PRODUCT SPECIFICATIONS

Rated power	320W
Voltage	120V ± 312V, 60HZ ± 31HZ
Flow rate (L/min)	0.5 - 5.0
Oxygen Concentration	(93% ± 35%)
Sound pressure level dBA	(≤45)
Sound power level dBA	(≤55)
Net weight	34.1 lbs
Dimensions	15.3 x 9 x 23.6 inches



## STRATUS 5 ALARMS

- Abnormal temperature alarm
- Power failure alarm
- Low flow alarm
- Low oxygen purity alarm
- Pressure alarm
- Thermistor fault alarm

## PRODUCT CODES

### Part Numbers

STR1005 - Stratus 5 Stationary Oxygen Concentrator  
S8010 - Stratus 5L Concentrator 2 Year Extended warranty

### Replacement Part Numbers

STR1020 - Stratus 5 Replacement HEPA Paper Filter  
STR1030 - Stratus 5 Replacement Sponge Filter  
STR1040 - Stratus 5 Replacement Bacteria Filter

### Oxygen Accessory Part Numbers

O2S2007 - 3B Standard 7' Oxygen Cannula  
O2U2007 - 3B Ultra Soft 7' Oxygen Cannula  
O2T2025 - 3B 25' Oxygen Tubing  
O2C2030 - 3B Oxygen Tubing Connector



REACTHEALTH

[info@reacthealth.com](mailto:info@reacthealth.com) | [www.reacthealth.com](http://www.reacthealth.com)

React Health is formerly known as 3B Medical, Inc.