React Health PAP Masks: Offering the right fit, the first time

SIESTA FULL FACE MASK



REFIT RATE AS EVIDENCED BY SLEEPGLAD AGGREGATE DATA



SIESTA NASAL MASK



REFIT RATE AS EVIDENCED BY SLEEPGLAD AGGREGATE DATA


RIO II NASAL PILLOWS MASK



REFIT RATE AS EVIDENCED BY
SLEEPGLAD
AGGREGATE DATA

2020- April 2021







Our Family Of CPAP Mask Interfaces



The Siesta

Full Face Mask

- No forehead brace allows for increased line of sight - even with your glasses on
- Self-adjusting ergonomic cushion design allows for superior seal
- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- 360-degree dual-swivel elbow for active sleepers and allows for sleep in any position
- Air vent diffuser for a quiet nights sleep
- Anti-phyxia valve safety feature opens for fresh air in case of machine shut off





The Siesta

Nasal Mask

- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow for simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- 360-degree dual-swivel elbow for active sleepers and allows for any sleep position



The Rio[™]II Nasal Pillows Mask

- A lightweight experience- Just 2.3 ounces: One of the lightest masks on the market today
- Dual wall pillows with trampoline bounce for secure fit
- Split back strap allows for a better fit, while
 Velcro fasteners create added layer of adjustment
- Unique rotating ball and socket swivel for active sleeper

