

INTRODUCING

Stratus 5

The 5 LPM Stationary Oxygen Concentrator from React Health

- Slim line design for a smaller footprint in the home
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet operation compressor is specially placed, and padded to reduce noise
- Simple, intuitive operation





Stratus 5L Designed with the Patient in Mind

STRATUS 5L

Designed with your business in mind

Simple to service
Easy access to internal parts

3-year standard warranty



PRODUCT SPECIFICATIONS

Rated power	320W
Voltage	120V ± 312V, 60HZ ± 31HZ
Flow rate (L/min)	0.5 - 5.0
Oxygen Concentration	(93% ± 35%)
Sound pressure level dBA	(≤45)
Sound power level dBA	(≤55)
Net weight	34.1 lbs
Dimensions	15.3 x 9 x 23.6 inches



STRATUS 5 ALARMS

- Abnormal temperature alarm
- Low oxygen purity alarm
- Power failure alarm
- Pressure alarm

- Low flow alarm
- Thermistor fault alarm

PRODUCT CODES

Part Numbers

STR1005 - Stratus 5 Stationary Oxygen Concentrator

Replacement Part Numbers

STR1020 - Stratus 5 Replacement HEPA Paper Filter

STR1030 - Stratus 5 Replacement Sponge Filter

STR1040 - Stratus 5 Replacement Bacteria Filter

Oxygen Accessory Part Numbers

O2S2007 - 3B Standard 7' Oxygen Cannula

O2U2007 - 3B Ultra Soft 7' Oxygen Cannula

O2T2025 - 3B 25' Oxygen Tubing

O2C2030 - 3B Oxygen Tubing Connector



info@reacthealth.com | www.reacthealth.com

React Health is formerly known as 3B Medical, Inc.